## **Mermaid Premium Salt**

## **Product Description**

Mermaid Premium Salt is a refined, dried and naturally evaporated sea salt (sodium chloride) suitable for human consumption complying with the Australian New Zealand Food Standards Code 2.10.2. This product is a translucent to white solid, with a texture of medium fine granules and a mean particle size of approx. 0.45 mm. It is produced in Australia by natural solar evaporation of seawater, harvested, washed, dried, sized and packed in accordance with good manufacturing practice, under a quality system that complies with ISO 9001.

## **Storage Conditions**

Product is shelf life stable. Long-term storage does not adversely affect salt except for caking or lumping as salt absorbs/expels moisture from/to the atmosphere. Fine grain salts are particularly susceptible to caking. As a guide this product should be used within 3 - 6 months from manufacture date to avoid significant caking problems, however customers should assess their own individual needs for ordering frequency, stock rotation, stock levels and local conditions. To avoid significant caking, adjust ordering volume and frequency.

Chemical Analysis	Specification	
Purity (%NaC1 min dry basis)	99.4	
Moisture (% max)	0.20	
Insolubles (% max)	0.03	
Magnesium (mg/kg max)	500	
Calcium (mg/kg max)	1000	
Sulphate (mg/kg max)	2500	
Iron (mg/kg max)	10	

Grainsize (Cumulative % Retained)					
	Min	Max			
0.85 mm	0%	1%			
1.71 mm	0%	25%			
0.212 mm	80%	100%			
0.150 mm	95%	100%			
Through 0.150mm not greater than 5%					

e (Cumulative % Retained)		etained)	Palletising - Standard
	Min	Max	Available in variety of bag sizes and pallet configurations
	0%	1%	Contains No Additives
	0%	25%	
	80%	100%	
	95%	100%	
0.150mm not greater than 5%		er than 5%	

## Protection and Labelling

Packed and Sealed in Woven Polypropylene Bags (with polyethylene liner)

Bonded Internal Laminate

Pallet Slip Sheet and Stretch Wrapped

Traceability: Date of Manufacture or Batch Number which includes DOM in format YYMMDD-BCH and then an automated sequential number ie YYMMDD-BCH-123456; printed on the side of individual bags as well as on pallet label

Method of Analysis					
Purity (%NaC1 min dry basis)	In house method (Ref. ASTM International Standard E534 Standard Test Methods for Chemical Analysis of Sodium Chloride)				
Moisture (% max)	In house method (Ref. ASTM International Standard E534 Standard Test Methods for Chemical Analysis of Sodium Chloride)				
Insolubles (% max)	In house method (Ref. ASTM International Standard E534 "Water Insolubles" 100g used, results reported to 0.001%)				
Calcium & Magnesium	Reference: ASTM E 534 "Calcium Magnesium" Used" Japan Tobacco & Salt: Methods for Salt Analysis Item 5.1 "Calcium & Magnesium" (atomic absorption)				
Calcium (mg/kg max)	Reference: ASTM E 534 "Calcium Magnesium" Used" Japan Tobacco & Salt: Methods for Salt Analysis Item 5.1 "Calcium & Magnesium" (atomic absorption)				
Sulphate (mg/kg max)	Reference "AS 2093.1977 Appendix M" Used: turbidimetric analysis on UV-vis spectrophotometer				
Iron (mg/kg max)	In house method (Ref. AACC) 40-41.03 Iron - Spectrophotometric Method				
Grainsize	ISO 2591-1:1988 (R2017) Test Sieving - Methods using test sieves of woven wire cloth and perforated metal plate				

Level 6, 565 Bourke Street, Melbourne VIC 3000 Ph: 1800 032 046 sales@cheethamsalt.com.au





Issued: 22 October 2024